



# Winkonlaser

Bringing Beauty and Love by Technology

Add: D2510, Victoria Building, Beijing, China

WhatsApp: +8618533632296

Email: [info@winkonlaser.com](mailto:info@winkonlaser.com)


Web: [www.winkonlaser.com](http://www.winkonlaser.com)/[www.beautemax.com](http://www.beautemax.com)

## Winkonlaser®

Bringing Beauty and Love by Technology

### RF RENASLIM® FAT

- Muscle Building
- Fat Reduction
- Skin Tightening





## Introducing **RENASLIM**<sup>®</sup>

**RENASLIM**<sup>®</sup> Fat Dissolving Machine, a state-of-the-art, radio frequency energy-based therapy with real-time temperature control. Clinically proven to permanently eliminate stubborn fat cells in areas such as abdomen, flanks, arms, bra straps, legs, double chin and knees. In addition to reducing fat, RF energy can also promote collagen regeneration, improve skin sagging and wrinkles, improve the overall contour of the face and jaw line, and make the face and body skin firm and smooth, reducing fat, shaping and firming.



The  
**RENASLIM**<sup>®</sup>  
Designed for  
you.

### ✓ Safe

Non-invasive and non-ablative. It is easy to use without an operator, and it is simple and safe.

### ✓ Versatile

Special hand-held handle, suitable for smaller and more sensitive parts of the body, such as side breasts, double chin, face. Maintain the subcutaneous temperature of the treated area evenly throughout the procedure.

### ✓ Easy to use

Simultaneous treatment of multiple areas. 6 (flat fixed) hands-free handles that can cover 300cm<sup>2</sup> at the same time on the abdomen and both sides.

### ✓ Smart Monitoring Software

Equipped with a new large sized 12 inch Android screen, with pre-set parameters, practitioners can enjoy a more intuitive, user friendly work experience.









# Synergetic Technologies. Superior Results.

The most special 6 flat fixed treatment handles and 1 handle sliding treatment handles break through the traditional radio frequency mobile operation, bringing a safe, efficient, fast and most comfortable experience, and can be easily used without an operator



Use 2MHZ mono polar radio frequency to deliver even energy to the fat layer and then maintain a comfortable skin temperature. Allows you to treat multiple areas of the body simultaneously for maximum results and patient satisfaction.



Real-Time Temperature Control: Continuously monitors skin temperature and automatically adjusts energy delivery to comfortably heat and maintain skin to 45° C .surface temperature cooler than deep fat temperatures.



Fat heating selectivity: Optimized energy delivery and heating enhance selective apoptosis within subcutaneous adipose tissue.



On average, 24-27% of fat cells are irreversibly damaged.




Over a 12-week period, damaged fat cells were slowly broken down and eliminated from the treated subcutaneous adipose tissue.



About 24 percent of the fat cells were irreversibly damaged 12 weeks after treatment, causing them to be excreted by the body.




# The **RENASLIM**<sup>®</sup>'s KEY FEATURES



Body contouring handle \*6  
Flat Fixation Therapy 50cm<sup>2</sup>/2MHz



Swan neck shaped  
handle \*1  
Double Chin Treatment  
16cm<sup>2</sup>



Sliding Shaping Handle \*1  
40cm<sup>2</sup>/2MHz



**Gift**  
Tie the skin better  
- Negative plate 4 pieces  
- Conductive stickers 20 pieces



Emergency handle \*1  
16cm<sup>2</sup>



# Clinical Case

