Bringing Beauty and Love by Technology

₩inkonlaser

Add: D2510, Victoria Building, Beijing, China

WhatsApp: +8618533632296

Email: info@winkonlaser.com

Web: www.winkonlaser.com/www.beautemax.com

5 IN 1

MULTIFUNCTION

- IPL
- SIFHR
- Nd:yag Laser
- Diode Laser
- RF



Introducing ML600

This is 5 in 1 Multifunctional beauty device. It combines 5 technologies, E-light/IPL & Speed Hair Removal & Nd:yag laser & 808nm Diode Laser & RF. It has 4 handle pieces: Speed Hair Removal handle and 808nm Diode Laser and Nd:yag laser handle and RF handle.



IPL/SHR utilizes the principle of absorption of light to works on hair follicle. The light will be absorbed by Melanin of hair, after the hair follicle be destroyed. Then it will fall off automatically and be removed permanently.Speed Hair Removal also could promote recombination of collagen on the dermis, make the collagen active, and reach function of skin rejuvenation and wrinkle removal.

The ML500 **System Highlights**



Nd:yag laser erupt high energy in an instant, after absorbing the laser energy, the pigment in the tattoo broke instantly, then becomes tiny pigment group, finally get out from body with metabolism. Meanwhile, it also can do black doll skin rejuvenation treatment, using activated carbons unique characteristics , penetrate into pores, absorbing grease and dirt. Laser light stimulate the regeneration of dermal collagen and the shrink of collagen fibers, to reach a whitening, pigment lightening, pores shrinking and blackheads removal effect, making skin recover young.



808nm Diode Laser damage the chromophores by selectively heating them while leaving surrounding tissue unharmed.

• Permanent hair reduction • For all skin types and skin textures • Including skin rejuvenation treatment



RF can promote subcutaneous collagen contraction, while cooling the skin surface, the dermis layer is heated and the epidermis maintains normal temperature, this time will produce two kinds of reactions: First, the skin dermis thickens, wrinkles follow Lightening or disappearing; second, the remodeling of the subcutaneous collagen, resulting in new collagen, the skin becomes more compact after a treatment.



- Pigmented Lesions
 Skin Remodeling
 Tattoo Removal
 - Acne Removal Carbon Peeling Hair Removal
 - Vascular Therapy



Synergetic Technologies. **Superior Results**.

- Alarm protection system for water flow and water température.
- Ensure a safe treatment.
- Perfect cooling system: Water, Air.Semiconductor, Display program.
- Easy operation, quick return.
- No down time, no pain.


Wrinkle Removal

SHR/IPL Handle 🖌


• UK imported Xe lamp, more than 300000 flashes • Pure saphire crystal, up to 10 times flashes per second














530nm-1200nm Pigment treatment



640nm-1200nm Hair removal




SHR/IPL Handle 🖌





Before



















Nd:Yag Laser Handle 🍾



Before



After



Before



After




Before



Before



After



Before



After



Before









Diode Laser Handle

HSPC Contact Cooling

HSPC is an advanced technology that continuously cools skin and facilitates range-controlled temperature throughout the entire treatment. HSPC minimizes epidermal risks while maintaining heat within the dermis where hair follicles are treated. Temperature remains controlled and constant which contributes to the diode' s efficiency. Patients can now enjoy a completely pain-free, cool and comfortable experience throughout the session.

minim

Lifespan	50,000,000 sh
Spot Size	12mm*12mm
Laser Power	Up 1600W



USA Coherent imported laser , high durability.

nots



Diode Laser Handle



Before



After



Before



...

After



Before



Before



After



Before



After



Before











Eyes around wrinkle removal, forehead wrinkle




Face and neck wrinkle removal & lifting



Body fat reduction, Skin tightening , Body slimming







Before



After



Before



After




Before




Before



After



Before



After



Before



After

