



Winkonlaser
Bringing Beauty and Love by Technology

Add: D2510, Victoria Building, Beijing, China

WhatsApp: +8618533632296

Email: info@winkonlaser.com

Web: www.winkonlaser.com/www.beautemax.com

Winkonlaser®

Bringing Beauty and Love by Technology


HIEMT+EMS+RF

Three therapies in a single treatment.

BUILD MUSCLE

SCULPT YOUR BODY

SKIN TIGHTENING



Workout Smarter Not Harder

Non-essential fat provides greater resistance to athletic motion thereby forcing one to increase the muscle force of contractions per given workout. This takes time, energy and a lot of patience, until one discovers Renasculpt.



High-Intensity Focused Electro-Magnetic Technology

HIEMT technology, the full name is High-Intensity Focused Electro-Magnetic, it induce short bursts of powerful muscle contractions that lead to increased muscle density, decreased volume, better definition, and improved tone. Also the energy will eliminates the fat cells, so not only for muscle building, but also for fat reduction. Clinical research shows that after a course treatment, it can increase 16% muscle and get 19% fat reduction.



Electro Muscle Stimulation

EMS is a type of electrotherapy stimulates a muscle contraction using electrical impulses, this stimulation creates muscle contractions that can be quick and frequent, fast with long pauses, so get good results on muscle building as well.



Radio Frequency

RF is a nonsurgical method of tightening your skin. The procedure involves using energy waves to heat the deep layer of your skin known as your dermis. This heat stimulates the production of collagen so that get skin firming and tightening results.



RF+EMS

HIEMT

Combines 3 energies for simultaneous fat elimination and muscle building in a 30-minute session.

Multiple working modes, Provide diverse Therapeutic experience
Two HI-EMT + EMS&RF working modes

- ① Combined mode, three technologies work at the same time.
- ② Alternate mode: HI-FEM first, then EMS&RF, cycle work



Before



After



Big Handle × 2

You begin to feel tangible results right after the treatment. Positive results are usually reported two to four weeks after the last session and continue to improve for several weeks following the treatments.



Small Handle × 2

ReNASCUPT treatment does not require any pre-treatment preparations. You can follow your regular lifestyle right after the treatment.



Pelvic Chair Cushion

- 1- The active regeneration of collagen in pelvic floor tissues, strengthens pelvic floor muscles, and lifts buttocks
- 2- Improve privacy issues, such as urinary incontinence, penetration, privacy issues
- 3- Achieve vagina tightening.




Before



After



Before



After



Product Name	FE80
Technology	HIEMT+RF+EMS
Handle Piece	2 Big Size & 2 Small Size & Pelvic Chair Cushion
Magnetic Vibration Intensity	15 Tesla
RF Frequency	20MHZ
HIFEM Frequency	1-135HZ
Treatment Mode	Combined Mode& Alternate Mode
Input Voltage	AC110V-220V
Output Power	5000W
Pulse Duration	300μs
Fuse	20A
Cooling System	Air Cooling
Size of Flight Shipping Case	69*57*112cm
Gross /Net Weight	90/71KG
Working Mode	Automatic Mode and Manual Mode